

# SPECTROLINE®

## OFK-450 LED Forensic Inspection Kit

Features the super-powerful, portable and rechargeable **OPTIMAX™ 450** blue light LED flashlight. Ideal for crime scene investigation, gathering potential evidence and laboratory work.

Enhances the detection, inspection and photography of:

- ✓ Bodily fluids
- ✓ Biological stains and latent fingerprints
- ✓ Bruises, bite marks, pattern wounds, hair and fibers



### FEATURES:

- Super-high intensity, 450 nm blue light LED
- 50,000-hour LED life
- Instant-on operation; lamp reaches full intensity immediately
- Corrosion-resistant, black-anodized aluminum body stands up to heavy use in the field or laboratory
- Detachable diffusing filter ensures uniformity of light source. Ideal for forensic photography.
- Powered by rechargeable NiMH battery. Provides 90 minutes of continuous inspection between charges.



**The OFK-450 Forensic Inspection Kit** comes complete with the OPTIMAX™ 450 blue light LED inspection flashlight, detachable diffusing filter, smart AC and DC chargers, orange contrasting spectacles, belt holster and padded carrying case. AC charger available in 120V, 230V, 240V or 100V versions.

**SPECTRONICS CORPORATION**  
956 Brush Hollow Road, P.O. Box 483  
Westbury, New York 11590  
800-274-8888 • 516-333-4840  
Fax: 800-491-6868 • 516-333-4859  
[www.spectroline.com](http://www.spectroline.com)

DISTRIBUTED BY:

ISO 9001:2000 CERTIFIED COMPANY

3/09 A09026  
PRINTED IN U.S.A.